

Medium Steam Sterilizers



Model Azteca AC Series	AC-450	AC-470	AC-575	AC-5100
Chamber dimensions, Inner Dia. x Depth	400 x 670 mm	400 x 750 mm	500 x 700 mm	500 x 1000 mm
Chamber volume	75 l	90 l	150 l	200 l
External dimensions W x H x D	750 x 630 x 970 mm		850 x 730 x 1050 mm	850 x 730 x 1300 mm
Approximate Weight	160 kg	170 kg	200 kg	240 kg
Power supply*	3 ph. 380-400 V, 50/60 Hz			
Power with steam generator	7 kW	9.9 kW	12 kW	
Available sterilization cycles	Flash 134 °C / Unwrapped 134 °C / Wrapped 134 °C / Prion 134 °C / Porous 121 °C			
Available test cycles	Bowie & Dick test / Vacuum test			
Touch-screen display	LCD 5.7" Color Graphic			
Cabinet	Painted quality steel		Stainless steel	

*Adjustable to different voltage system